



RamanChada

Cheat Sheet

Loading & displaying data

R = **RamanChada**('C:\file.spc')

Spectrum object info: R

RamanChada with 3526 points generated Wed Sep 15 ...

Get metadata as dict: R.meta

Add metadata (dict):

R.**add_metadata**({'Power[mW]':5})

Get processing log as list:

R.**log** or R.**show_log**()

Undo last processing step: R.**rewind**(-1)

Revert to original data: R.**rewind**(0)

Save .cha file with all changes:

R.**commit**('commit message')

Load .cha file:

R = **RamanChada**('C:\file.cha')

Plot spectrum: R.**plot**()

Get the raw data from the same file:

S = **RamanChada**(
R.file_path, raw=True)

Cosmic rays


Fit x ray model

R.**fit_xrays**()

Plot x ray model

R.**plot_xrays**()

Remove x rays: R.**remove_xrays**()



Calibration

Interpolate to x axis of reference spectrum:

R.**interpolate_x**(reference_spectrum)

Calibrate x with existing RamanCalibration cal:


R.**calibrate**(cal)

Show x calibration curve:

cal.**show**()

Get calibration time / date:

cal.time



Generate x calibration using reference spec:

cal = R.**make_x_calibration**(ref)

...and consider only peaks in interval:

cal = R.**make_x_calibration**(
ref,1700,3200)

Generate x calibration using peak positions list:

cal = R.**make_x_calibration**(
[202.12,451.76,...,1809.28])

Calibrate y with existing RamanCalibration:

R.**calibrate_y**(y_cal)

Generate y calibration using reference spectrum:

y_cal = R.**make_y_calibration**(ref)

Save RamanCalibration to disk:

cal.**save**('C:\cal_filename.chacal')

Load RamanCalibration from disk:

cal = **read_x_calibration**(
'C:\cal_filename.chacal')

Baseline separation

Fit baseline using SNIP method


R.**fit_baseline**(method='snip')

Plot baseline model

R.**plot_baseline**()

Fit and remove baseline

R.**remove_baseline**(method='snip')



Pre-processing

Smooth spectrum using Savitzky-Golay filter:

R.**smooth**('sg',window=11,order=3)

Normalize spectrum using vector norm:

R.**normalize**('vector')

Area normalization using only an interval:

R.**normalize**('area',500,1250)

Crop spectrum on the x axis to 500-1250 cm⁻¹:

R.**x_crop**(500,1250)

Add spectrum S to R: R.**math**(S, '+')

Peaks search & fitting

Find peaks with prominence>0.2 without fitting:

R.**peaks**(prominence=0.2,fit=False)

DataFrame with detected/fitted peaks: R.bands

Fit peaks within 2xFWHM w/ Voigt profile & plot:


R.**peaks**(fitmethod='voigt',
interval_width=2,show=True)

Find peaks using wavelets, fit & sort by position:

R.**peaks**(cwt=True,sort_by='position')

*Plot spectrum with
peak positions*

R.**show_bands**()



Batch processing

List of RamanChada objects from path list files:

SL = [**RamanChada**(f) for f in files]

Normalize all spectra by standard normal variate:

[s.**normalize**() for s in SL]

Apply Wiener filter to all and then save to disk:

[s.**smooth**('wiener',7) for s in SL]

[s.**commit**('smoothed') for s in SL]